



CLINICAL NOTES

ADDING CLINICAL NOTES

Clinical notes may now be added. They will appear at the end of the last page on the patient's flowsheet report.

To create new notes:

1. Click the desired patient cell from the main Monitor screen. A detail view will appear. Select the Notes tab located on the right hand side of the screen.
2. Select the New button to enter new clinical notes to the patient's detail cell. These will appear at the end of the flowsheet report.
3. Enter clinical notes and your initials. Click Save when finished.
4. Once your notes are saved, you may view other notes for the patient by toggling between them with the Next and previous buttons.



The Close button will return you to the patient's detail screen. Select the Monitor button to return to the main screen.

CREATING CUSTOM FIELDS

Users who want custom fields on their reports should contact their Bernoulli system administrator. From the Bernoulli control panel, users can create custom parameters and also designate specific colors for these parameters. Please refer to the Bernoulli Administrators manual for using the control panel and customizing flowsheet reports.



Cardiopulmonary Corporation
200 Cascade Blvd.
Milford, CT 06460
(203) 877-1999; Fax: (203) 877-3401
(800) 337-9936
www.cardiopulmonarycorp.com

CAUTION: UNITED STATES FEDERAL LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.
NOTICE: United States Patents 6,839,753, and International Equivalents. Additional U.S. and International Patents pending.
Bernoulli® is a US registered trademark of Cardiopulmonary Corp.
© 2006 Cardiopulmonary Corp. All Rights Reserved.

BERNOULLI® FLOWSHEETS


QUICK START GUIDE



PART NO. 2741 | NOVEMBER 2006



FLWSHEET REPORTING



1. Select Report from the Main Screen. (shown above)

2. Select Flowsheets from the Report main menu.(fig.1)




fig.1

3. Choose desired patient, admitted or recently discharged, by selecting the box next their name and room number. Then select Continue with selected patients when finished. (fig.2)




fig.2

4. Select the information to be shown on the flowsheet report from the Flowsheet options screen (fig.3). This screen allows you to select the types of items displayed on the snapshots list. Each column in the report will represent the values at the time a Snapshot was taken.

- ▶ Select **Scheduled** to include scheduled snapshots
- ▶ Select **Manual** to include manual snapshots.
- ▶ Select **Setting Change** to include snapshots triggered by setting changes.
- ▶ Select **Alarm** to include snapshots triggered by alarms.




fig.3

When finished select Use these options.

5. The Time Range Selection Screen will appear. To choose a time range, select one of the "most recent" time range buttons offered or use the + and - buttons to set a start and stop time. Click Use this Time Range when finished. (fig.4)





fig.4



6. The Flowsheet Report will appear. (shown above)

Using the Controls:

- ▶ Click the Zoom Out and Zoom In buttons to minimize or increase the Flowsheet view.
- ▶ Use Next Page and Prev Page to toggle between flowsheet page views.
- ▶ Clicking First Page will bring you back to the beginning of the flowsheet.
- ▶ Select Last Page to view the end of the flowsheet document. This is where clinical notes can be viewed. To add notes see reverse.
- ▶ Select the Pages side tab to see a preview of what the flowsheet will look like when printed.

UPDATED FEATURES

A. Color fields have been added for alarm Violations. Measurements in alarm will be shown with a red cell and black text. Measurements that are 10% away from an alarm limit and not in alarm will be shown in a yellow cell with black text. Since the flowsheet may be printed in black and white, alarm functions will be printed in bold with asterisks signifying the corresponding alarm measurement.

B. Custom headers, such as a hospital logo, may be added to the flowsheet for each customer.

C. Bold indications on the flowsheet have now been limited to setting changes in order to catch the eye of the reader. The previous system used bolding on all the constantly changing patient and device values. (Fig.3)